

Table 1 Comparison between enzyme-linked immunosorbent assay (ELISA) IgM and IgG-avidity results in the detection of past or recent toxoplasmosis according to polymerase chain reaction (PCR) results in samples from pregnant women at risk of toxoplasmosis (n = 25 samples)

IgG avidity	ELISA results		Total No.	PCR results	
	IgM status			DNA +ve (n = 16) No.	DNA -ve (n = 9) No.
Low	-ve		0	0	0
	+ve		15	15	0
High	-ve		7	0	7
	+ve		3	1	2
Total			25	16	9

+ve = positive; -ve = negative.